




Medical Blue Sensitive X-Ray Film



Product Features

- This film is used for the examination of human body disease and abnormal changes of viscera, bones and soft tissues, providing reliable information for clinical diagnosis.
- The film has sensitive emulsion layers on both sides of bluish polyester base, and must be used with the blue emitting intensifying screen.
- It features high sensitive speed, moderate contrast, abundant gradation, good definition and low fog.

To learn more, email us: gracewtj@gmail.com or call (0086) 755-8146 3292

Handling

- Can be matched with blue intensifying screen or the screen of the other same system which matches the blue sensitive x-ray.
- The film should be handled carefully by the edges to avoid physical strains.
- The best conditions of usage should be: temperature:18-24°C, humidity :50%RH-65%RH.
- Hands must be clean, dry and free of lotions ,etc.

Film Processing

- This film is designed for 90-seconds automated processing or manual processing.
- The former is recommended with any kinds of film processor, developing at 30-35°C in 23s.
- In manual processing ,the film should be developed in temperature 20°C for 5 minutes, fixed in 5-8 minutes, washed for 10-30 minutes and dried.

Packing and specifications

- **Export Packing:** paper box, 5boxes/ carton; 60 cartons/Pallet(12*15in)
- **Outer Packing Size:**

18*24cm	100 NIF	5boxes	295*180*235mm	0.01248 m3
24*30cm	100 NIF	5boxes	340*180*300mm	0.01836 m3
30*40cm	100 NIF	5boxes	440*180*365mm	0.02891 m3
35*35cm	100 NIF	5boxes	400*180*405mm	0.02916 m3
35*43cm	100 NIF	5boxes	480*180*405mm	0.03499 m3

Mark: for a 20" container, it can put about 4,200boxes 35*35cm.

- **Shelf Life:** 24months

To learn more, email us: gracewtj@gmail.com or call (0086) 755-8146 3292


Safelight

- Use a safelight filter (dark red) in a suitable safelight lamp equipped with a 15-watt bulb. Keep the film at least 4 feet (1.2 meters) from the safelight, and the total lighting time is in 10 minutes.

Storage

- Store in a cool, dry place (The suitable temperature is 10-20°C, the humidity is 30%RH-60%RH).
- Keep it away from dampness, high temperature, any bad effective gases, x-ray, gamma ray or any other penetrating rays.
- Stored vertically to avoid any harmful effect caused by pressure.

Characteristic Curve





Medical Green Sensitive X-Ray Film



Product Features

- This film is used for the examination of human body disease and abnormal changes of viscera, bones and soft tissues, providing reliable information for clinical diagnosis.
- The film is a general green sensitive medical x-ray film, which has green sensitive emulsion layer on both sides of bluish polyester bases, and must be used with the general green emitting intensifying screen.
- It features high sensitive speed, good definition, wide exposure latitude and good using compatibility.

To learn more, email us: gracewtj@gmail.com or call (0086) 755-8146 3292

Handling

- Can be matched with green intensifying screen or the screen of the other same system which matches the green sensitive x-ray.
- The film should be handled carefully by the edges to avoid physical strains.
- The best conditions of usage should be: temperature:18-24°C, humidity :50%RH-65%RH.
- Hands must be clean, dry and free of lotions ,etc.

Film Processing

- This film is designed for 90-seconds automated processing or manual processing.
- The former is recommended with any kinds of film processor, developing at 30-35°C in 23s.
- In manual processing ,the film should be developed in temperature 20°C for 5 minutes, fixed in 5-8 minutes, washed for 10-30 minutes and dried.

Packing and specifications

- **Export Packing:** paper box, 5boxes/ carton; 60 cartons/Pallet(12*15in)

- **Outer Packing Size:**

8*10in	100 NIF	5boxes	300*180*255mm	0.01377 m3
10*12in	100 NIF	5boxes	350*180*310mm	0.01953 m3
11*14in	100 NIF	5boxes	400*180*325mm	0.02340 m3
12*15in	100 NIF	5boxes	425*180*350mm	0.02678 m3
14*14in	100 NIF	5boxes	400*180*405mm	0.02916 m3
14*17in	100 NIF	5boxes	480*180*405mm	0.03499 m3

Mark: for a 20" container, it can put about 4,200boxes 35*35cm.

- **Shelf Life:** 24months

To learn more, email us: gracewtj@gmail.com or call (0086) 755-8146 3292

Safelight

- Use a safelight filter (dark red) in a suitable safelight lamp equipped with a 15-watt bulb. Keep the film at least 4 feet (1.2 meters) from the safelight, and the total lighting time is in 10 minutes.

Storage

- Store in a cool, dry place (The suitable temperature is 10-20°C, the humidity is 30%RH-60%RH).
- Keep it away from dampness, high temperature, any bad effective gases, x-ray, gamma ray or any other penetrating rays.
- Stored vertically to avoid any harmful effect caused by pressure.

Characteristic Curve

